Technical Document

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
<参考文献>

1. Ishihara K, Ueda T, Nakabayashi N. : Preparation of phospholipid polymers and their properties as polymer hydrogel membrane. Polym J. 22, 355-360,1990.
2. Ishihara K, Iwasaki Y, Ebihara S, Shindo Y, Nakabayashi N. : Photoinduced graft polymerization of 2-methacryloyloxyethyl phosphorylcholine on polyethylene membrane surface for obtaining blood cell adhesion resistance. Colloids Surf B Biointerfaces. 18, 325-335, 2000.
3. Moro T, Kawaguchi H, Ishihara K, Kyomoto M, Karita T, Ito H, Nakamura K, Takatori Y. : Wear resistance of artificial hip joints with poly(2-methacryloyloxyethyl phosphorylcholine) grafted polyethylene: comparisons with the effect of polyethylene cross-linking and ceramic femoral heads. Biomaterials. 30(16), 2995-3001, 2009.
4. Moro T, Takatori Y, Ishihara K, Nakamura K, Kawaguchi H. : 2006 Frank Stinchfield Award: grafting of biocompatible polymer for longevity of artificial hip joints. Clin Orthop Relat Res. 453, 58-63, 2006.
5. Kyomoto M, Moro T, Takatori Y, Kawaguchi H, Ishihara K. : Cartilage-mimicking, High-density Brush Structure Improves Wear Resistance of Crosslinked Polyethylene: A Pilot Study. Clin Orthop Relat Res. 469, 2327-2336, 2011.
6. Kyomoto M, Moro T, Konno T, Takadama H, Kawaguchi H, Takatori Y, Nakamura K, Yamawaki N, Ishihara K. : Effects of photo-induced graft polymerization of 2-methacryloyloxyethyl phosphorylcholine on physical properties of cross-linked polyethylene in artificial hip joints. J Mater Sci Mater Med. 18(9), 1809-15, 2007.
7. Moro T, Takatori Y, Ishihara K, Konno T, Takigawa Y, Matsushita T, Chung UI, Nakamura K, Kawaguchi H. : Surface grafting of artificial joints with a biocompatible polymer for preventing periprosthetic osteolysis. Nat Mater. 3(11), 829-36, 2004.

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